





PERIPHERAL INTERVENTION




Dredger[™]
Supporting Catheter

DREDGER™: THE TRUSTWORTHY CHOICE FOR DRILL AND DREDGE IN ARTERIOSCLEROSIS OBLITERANS




Dual Tapered Tip

- Excellent tip profile for easier lesion penetration




3 Radiopaque Markers + High Radiopaque Polymer Tip

- Precise catheter positioning
- Easy to estimate lesion length




Transparent Tube Design

- Visible of blood



Double Flat Wire Braided Structure

- Provides with strong support and kink resistance



45° 30° 17° Straight

4 Pre-shaped Curves

- Suitable for different vessel morphology
- Assist in guidewire control and lesion opening

ORDERING INFORMATION

Standard version*											
Guide wire (max)	Tip O.D (mm/Fr)	Proximal O.D (mm/Fr)	Distal O.D (mm/Fr)	Angiography Catheter (Fr)	Sheath (Fr)	Curve	Effective Length (cm)				
							65	90	135	150	170
0.014"	0.43/1.3	0.9/2.7	0.7/2.1	4	4	Straight	1681406500	1681409000	1681413500	1681415000	1681417000
						17°	1681406517	1681409017	1681413517	1681415017	1681417017
						30°	1681406530	1681409030	1681413530	1681415030	1681417030
						45°	1681406545	1681409045	1681413545	1681415045	1681417045
0.018"	0.54/1.6	0.95/2.9	0.8/2.4	4	4	Straight	1681806500	1681809000	1681813500	1681815000	1681817000
						17°	1681806517	1681809017	1681813517	1681815017	1681817017
						30°	1681806530	1681809030	1681813530	1681815030	1681817030
						45°	1681806545	1681809045	1681813545	1681815045	1681817045
0.035"	1.00/3.0	1.4/4.2	1.3/3.9	6	4	Straight	1683506500	1683509000	1683513500	1683515000	1683517000
						17°	1683506517	1683509017	1683513517	1683515017	1683517017
						30°	1683506530	1683509030	1683513530	1683515030	1683517030
						45°	1683506545	1683509045	1683513545	1683515045	1683517045

Expanded version**											
Guide wire (max)	Tip O.D (mm/Fr)	Proximal O.D (mm/Fr)	Distal O.D (mm/Fr)	Angiography Catheter (Fr)	Sheath (Fr)	Curve	Effective Length (cm)				
							65	90	135	150	170
0.018"	0.54/1.6	1.05/3.2	0.8/2.4	5	4	Straight	1681806500L	1681809000L	1681813500L	1681815000L	1681817000L
						17°	1681806517L	1681809017L	1681813517L	1681815017L	1681817017L
						30°	1681806530L	1681809030L	1681813530L	1681815030L	1681817030L
						45°	1681806545L	1681809045L	1681813545L	1681815045L	1681817045L

Frequently used
 Commonly used

- The 0.035" standard version can coaxially match with the 0.014"/0.018" standard version* seamlessly
- 5F angiography catheter can coaxially match with the 0.018" expanded version** seamlessly



@APT Medical

Address: Room 601, Building B, Tongfang Information Port, No.11 Langshan Road, Nanshan District, Shenzhen, P.R. China
 Email: international@aptmed.com Tel: +86 755-83591013
 PN: ENG-Supporting catheter-210X285X4P-20240313

